

1117 Low Alloy Steel Bar 41mm D x 150mm

Art. ID BS-03D
Unit bar
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | <0.001 | % | | |
| | Carbon (C) | [7440-44-0] | 0,18 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,02 | % | | |
| | Chromium (Cr) | [7440-47-3] | 0,18 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,27 | % | | |
| | Manganese (Mn) | [7439-96-5] | 1,15 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 0,04 | % | | |
| | Nitrogen (N) | [7727-37-9] | 0,011 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,11 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,025 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,1 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,28 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,017 | % | | |