

SV Mix, CLP method - Acid Surrogate Std 3/90 Statement of work (SOW) 1500µg/mL, Methanol

Art. ID RES-31003
Unit 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Fluorophenol	[367-12-4]	1500	µg/mL		
	Phenol D6	[13127-88-3]	1500	µg/mL		
	2-Chlorophenol-3,4,5,6-D4	[93951-73-6]	1500	µg/mL		
	2,4,6-Tribromophenol	[118-79-6]	1500	µg/mL		