

**SV Mix, CLP method - Internal Std 4000 µg/mL, CH<sub>2</sub>CL<sub>2</sub>**

Art. ID RES-31006  
Unit 1 mL ampule  
Deliverydetails De Minimis UN1593 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,4-Dichlorobenzene D4	[3855-82-1]	4000	µg/mL		
	Naphthalene D8	[1146-65-2]	4000	µg/mL		
	Acenaphthene D10	[15067-26-2]	4000	µg/mL		
	Phenanthrene D10	[1517-22-2]	4000	µg/mL		
	Chrysene D12	[1719-03-5]	4000	µg/mL		
	Perylene D12	[1520-96-3]	4000	µg/mL		