

SV Mix, 600 series method - 604 Calibration Std phenols 2000 µg/mL, methanol

Art. ID RES-31029
Unit 1 mL ampoule
Deliverydetails De Minimis UN1992 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phenol	[108-95-2]	2000	µg/mL		
	2-Chlorophenol	[95-57-8]	2000	µg/mL		
	2-Nitrophenol	[88-75-5]	2000	µg/mL		
	2,4-Dimethylphenol	[105-67-9]	2000	µg/mL		
	2,4-Dichlorophenol	[120-83-2]	2000	µg/mL		
	4-Chloro-3-methylphenol	[59-50-7]	2000	µg/mL		
	2,4,6-Trichlorophenol	[88-06-2]	2000	µg/mL		
	2,4-Dinitrophenol	[51-28-5]	2000	µg/mL		
	4-Nitrophenol	[100-02-7]	2000	µg/mL		
	DNOC	[534-52-1]	2000	µg/mL		
	Pentachlorophenol	[87-86-5]	2000	µg/mL		