

**SV Mix,8000 series method - 8040 Calibration Std #1 phenols mix
2000µg/mL,Isopropanol,1ml/ampule**

Art. ID RES-31088
Unit 1 mL
Deliverydetails De Minimis UN1219 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phenol	[108-95-2]				
	3-Methylphenol	[108-39-4]				
	2-Nitrophenol	[88-75-5]				
	2,4-Dichlorophenol	[120-83-2]				
	4-Chloro-3-methylphenol	[59-50-7]				
	2,4,6-Trichlorophenol	[88-06-2]				
	4-Nitrophenol	[100-02-7]				
	DNOC	[534-52-1]				
	Pentachlorophenol	[87-86-5]				