

SV Mix,8000 series method - 8040 Surrogate Std 2000µg/mL,Isopropanol,1ml/ampule

Art. ID RES-31090
Unit 1 mL
Deliverydetails De Minimis UN1219 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------|------------|---------------------|------|--------|--------|
| | 2-Fluorophenol | [367-12-4] | | | | |
| | 2,4,6-Tribromophenol | [118-79-6] | | | | |