

**QuEChERS Performance - Standard C, 300µg/mL Acetonitrile/Acetic Acid (99.0:0.1), 1mL/ampule**

Art. ID RES-31155  
Unit 1 mL  
Deliverydetails De Minimis UN1648 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bifenthrin	[82657-04-3]	300	µg/mL		
	Captan	[133-06-2]	300	µg/mL		
	Carbaryl	[63-25-2]	300	µg/mL		
	Cyprodinil	[121552-61-2]	300	µg/mL		
	Deltamethrin	[52918-63-5]	300	µg/mL		
	Fenhexamid	[126833-17-8]	300	µg/mL		
	Fenpropathrin	[39515-41-8]	300	µg/mL		
	Folpet	[133-07-3]	300	µg/mL		
	Imazalil	[35554-44-0]	300	µg/mL		
	Iprodione	[36734-19-7]	300	µg/mL		
	Metalaxyl	[57837-19-1]	300	µg/mL		
	Methiocarb	[2032-65-7]	300	µg/mL		
	Myclobutanil	[88671-89-0]	300	µg/mL		
	cis-Permethrin	[61949-76-6]	300	µg/mL		
	trans-Permethrin	[61949-77-7]	300	µg/mL		
	Thiabendazole	[148-79-8]	300	µg/mL		
	Vinclozolin	[50471-44-8]	300	µg/mL		