

### **Solder, 40Sn-60Pb (powder)**

Art. ID                      NIST-127b  
Unit                         150 g  
Deliverydetails            No Dangerous Good /not restricted

#### **Description**

This Standard Reference Material (SRM®) is in the form of powder. Certified values /// Sample value(s) - please ask for current certificate.

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Tin (Sn)          | [7440-31-5] | 39,3                | %    |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 0,43                | %    |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 0,01                | %    |        |        |
|                  | Bismuth (Bi)      | [7440-69-9] | 0,06                | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 0,011               | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,012               | %    |        |        |
|                  | Silver (Ag)       | [7440-22-4] | 0,01                | %    |        |        |