

Silicon Bronze (chip form)

Art. ID NIST-158a
Unit 150 g
Deliverydetails No Dangerous Good /not restricted

Description

Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	90,93	%		
	Silicon (Si)	[7440-21-3]	3,026	%		
	Zinc (Zn)	[7440-66-6]	2,076	%		
	Iron (Fe)	[7439-89-6]	1,228	%		
	Manganese (Mn)	[7439-96-5]	1,112	%		
	Tin (Sn)	[7440-31-5]	0,96	%		
	Aluminium (Al)	[7429-90-5]	0,4580	%		
	Lead (Pb)	[7439-92-1]	0,0973	%		
	Phosphorus (P)	[7723-14-0]	0,0263	%		