

## **Stainless Steel, Low-Carbon (AISI 316L) (chips)**

Art. ID	NIST-166c
Unit	100 g
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM®) is intended primarily for use in validation of chemical and instrumental methods of analysis. A unit of NIST-166c consists of a bottle containing approximately 100 g of chips. Certified value /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,00781	%		