

416 Se Stainless Steel Disc 50mm D x ~7mm XRF CUT

Art. ID BS-151-7MM
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | ~0,002 | % | | |
| | Carbon (C) | [7440-44-0] | 0,09 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,018 | % | | |
| | Chromium (Cr) | [7440-47-3] | 13,19 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,11 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,41 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 0,088 | % | | |
| | Nitrogen (N) | [7727-37-9] | 0,022 | % | | |
| | Niobium (Nb) | [7440-03-1] | 0,005 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,24 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,021 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,018 | % | | |
| | Selenium (Se) | [7782-49-2] | 0,328 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,65 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,005 | % | | |
| | Titanium (Ti) | [7440-32-6] | (<0,003) | % | | |
| | Vanadium (V) | [7440-62-2] | 0,046 | % | | |
| | Tungsten (W) | [7440-33-7] | 0,01 | % | | |