

Beryllium (Be), 1000 µg/mL in 2% H₂SO₄

Art. ID AP-CZ90071S1
Unit 100 mL/HDPE
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|-----------------|
| | Beryllium (Be) | [7440-41-7] | 1000 | µg/mL | | BeO: 99,99 % |