

Bismuth (Bi), 10000 µg/mL in 5% HNO3

Art. ID AP-CZ900810N1
Unit 100 mL/HDPE
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|-----------------|
| | Bismuth (Bi) | [7440-69-9] | 10000 | µg/mL | | Bi: 99,999 % |