

Copper (Cu), 1000 µg/mL in 2% HCl

Art. ID AP-CZ90151C5
Unit 500 mL/HDPE
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|-----------------|
| | Copper (Cu) | [7440-50-8] | 1000 | µg/mL | | Cu: 99,999 % |