

Iron (Fe), 1000 µg/mL in 2% HNO₃

Art. ID AP-CZ90191N5
Unit 500 mL/HDPE
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|-----------------|
| | Iron (Fe) | [7439-89-6] | 1000 | µg/mL | | Fe: 99,999 % |