

Stainless steel, Cr-Ni-Se (SAE 303Se)

Art. ID	NIST-339
Unit	150 g
Deliverydetails	No Dangerous Good /not restricted

Description

Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,052	%		
	Manganese (Mn)	[7439-96-5]	0,738	%		
	Phosphorus (P)	[7723-14-0]	0,129	%		
	Sulfur (S)	[7704-34-9]	0,013	%		
	Silicon (Si)	[7440-21-3]	0,654	%		
	Copper (Cu)	[7440-50-8]	0,199	%		
	Nickel (Ni)	[7440-02-0]	8,89	%		
	Chromium (Cr)	[7440-47-3]	17,42	%		
	Vanadium (V)	[7440-62-2]	0,058	%		
	Molybdenum (Mo)	[7439-98-7]	0,248	%		
	Cobalt (Co)	[7440-48-4]	0,096	%		
	Selenium (Se)	[7782-49-2]	0,247	%		