

## Ductile cast iron

Art. ID                      NIST-341  
Unit                         150 g  
Deliverydetails            No Dangerous Good /not restricted

### Description

Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	C total		1,81	%		
	C graphite		1,23	%		
	Manganese (Mn)	[7439-96-5]	0,024	%		
	Phosphorus (P)	[7723-14-0]	0,024	%		
	Sulfur (S)	[7704-34-9]	0,007	%		
	Silicon (Si)	[7440-21-3]	2,44	%		
	Copper (Cu)	[7440-50-8]	0,152	%		
	Nickel (Ni)	[7440-02-0]	20,32	%		
	Chromium (Cr)	[7440-47-3]	1,98	%		
	Vanadium (V)	[7440-62-2]	0,012	%		
	Molybdenum (Mo)	[7439-98-7]	0,01	%		
	Titanium (Ti)	[7440-32-6]	0,018	%		
	Magnesium (Mg)	[7439-95-4]	0,068	%		