

### **430F Stainless Steel Disc 38mm D x ~7mm XRF CUT**

Art. ID                      BS-154-7MM  
Unit                         disc  
Deliverydetails            No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Aluminium (Al)    | [7429-90-5] | ~0,002              | %    |        |        |
|                  | Carbon (C)        | [7440-44-0] | 0,03                | %    |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 0,019               | %    |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 17,58               | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 0,063               | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,4                 | %    |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7] | 0,31                | %    |        |        |
|                  | Nitrogen (N)      | [7727-37-9] | 0,039               | %    |        |        |
|                  | Niobium (Nb)      | [7440-03-1] | 0,005               | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,25                | %    |        |        |
|                  | Oxygen (O)        | [7782-44-7] | 0,008               | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,027               | %    |        |        |
|                  | Sulfur (S)        | [7704-34-9] | 0,302               | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 1,26                | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | ~0,005              | %    |        |        |
|                  | Vanadium (V)      | [7440-62-2] | 0,046               | %    |        |        |
|                  | Tungsten (W)      | [7440-33-7] | ~0,01               | %    |        |        |