

**Nickel (Ni), 1000 µg/mL in 2% HNO<sub>3</sub>**

Art. ID AP-CZ90381N5  
Unit 500 mL/HDPE  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit  | Method | Source          |
|------------------|-------------------|-------------|---------------------|-------|--------|-----------------|
|                  | Nickel (Ni)       | [7440-02-0] | 1000                | µg/mL |        | Ni: 99,999<br>% |