

Solva Tubing - 3 Tag, 95mm, Tag Colour: orange/green, ID: 0,38 mm, Wall: 0,91, OD: 2,2 mm

Art. ID	GE-0.38-OG3B-95-SF
Unit	pack of 12 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Tygon 4040 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Special tubing for hydrocarbons, petroleum products and distillates // Limitations: Not recommended for strong acids and alkalis, foodstuffs, beverages and medicines. Potential leaching of plasticizers // Advantages: Specially formulated to transport hydrocarbons, petroleum products and distillates, ideal for gasoline, kerosene heating oils, cutting liquids and coolants based on glycols, low gas permeability