

**Solva Tubing - 3 Tag, 95mm, Tag Colour: orange/yellow, ID: 0,51 mm, Wall: 0.86, OD: 2,23 mm**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | GE-0.51-OY3B-95-SF                |
| Unit            | pack of 12 tubes                  |
| Deliverydetails | No Dangerous Good /not restricted |

**Description**

All tubes made from Tygon 4040 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Special tubing for hydrocarbons, petroleum products and distillates // Limitations: Not recommended for strong acids and alkalis, foodstuffs, beverages and medicines. Potential leaching of plasticizers // Advantages: Specially formulated to transport hydrocarbons, petroleum products and distillates, ideal for gasoline, kerosene heating oils, cutting liquids and coolants based on glycols, low gas permeability