

Viton Tubing - 2 Tag, 127mm, Tag Colour: orange/white, ID: 0.64 mm, Wall: 0.90, OD: 2.44 mm

Art. ID	GE-VS0.64-OW
Unit	pack of 12 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Fluran HC F-5500A or equivalent // Physical Properties: fluoropolymerelastomer // Special Properties: Special tubing for concentrated acids and corrosive solvents // Limitations: Limited service-life // Advantages: High chemical resistance, low gas permeability, wide temperature range