

**Viton Tubing - 2 Tag, 95mm, Tag Colour: red/red, ID: 1.14 mm, Wall: 0.90, OD: 2.94 mm**

Art. ID	GE-VS1.14-RED-95
Unit	pack of 12 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Fluran HC F-5500A or equivalent // Physical Properties: fluoropolymerelastomer // Special Properties: Special tubing for concentrated acids and corrosive solvents // Limitations: Limited service-life // Advantages: High chemical resistance, low gas permeability, wide temperature range