

**Contour Flared End PVC Tubing - 2 Tag, 152mm, Tag Colour: orange/black, ID: 0.13 mm, Wall: 0.91, OD: 1.95 mm**

Art. ID	GE-0.13-OBLK-F
Unit	pack of 6 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Tygon R3603/R3607 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Inexpensive all-round tubing for general laboratory applications // Limitations: Potential leaching of plasticizers // Advantages: Transparent, resistant to almost all inorganic chemicals, does not age, low gas permeability, smooth inner surface