

Common Accessories - Contour Flared End PVC Tubing - 2 Tag, 152mm, Tag Colour: red/red, ID: 1,14 mm, Wall: 0,86, OD: 2,86 mm

Art. ID	GE-1.14-RED-F
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

Instruments: AGILENT (VARIAN) PERKIN ELMER SPECTRO LEEMAN JOBIN YVON /// All tubes made from Tygon R3603/R3607 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Inexpensive all-round tubing for general laboratory applications // Limitations: Potential leaching of plasticizers // Advantages: Transparent, resistant to almost all inorganic chemicals, does not age, low gas permeability, smooth inner surface