

Common Accessories - Contour Flared End PVC Tubing - 2 Tag, 152mm, Tag Colour: blue/blue, ID: 1.65 mm, Wall: 0.86, OD: 3.37 mm

Art. ID	GE-1.65-BLU-F
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Tygon R3603/R3607 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Inexpensive all-round tubing for general laboratory applications // Limitations: Potential leaching of plasticizers // Advantages: Transparent, resistant to almost all inorganic chemicals, does not age, low gas permeability, smooth inner surface /// Agilent (Varian)|MP-AES|4100/4200 MP-AES|Radial