

**Contour Flared End PVC Tubing - 2 Tag, 72mm, Tag Colour: orange/orange, ID: 0,89 mm, Wall: 0,91, OD: 2,71 mm**

Art. ID	GE-0.89-ORG-F-S
Unit	pack of 6 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Tygon R3603/R3607 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Inexpensive all-round tubing for general laboratory applications // Limitations: Potential leaching of plasticizers // Advantages: Transparent, resistant to almost all inorganic chemicals, does not age, low gas permeability, smooth inner surface