

Contour Flared End PVC Tubing - 3 Tag, 72mm, Tag Colour: purple/purple, ID: 2,06 mm, Wall: 0.86, OD: 3,78 mm

| | |
|-----------------|-----------------------------------|
| Art. ID | GE-2.06-PUR3B-F |
| Unit | pack of 6 tubes |
| Deliverydetails | No Dangerous Good /not restricted |

Description

All tubes made from Tygon R3603/R3607 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Inexpensive all-round tubing for general laboratory applications // Limitations: Potential leaching of plasticizers // Advantages: Transparent, resistant to almost all inorganic chemicals, does not age, low gas permeability, smooth inner surface