

**Contour Flared End PVC Tubing - 3 Tag, 72mm, Tag Colour: purple/black, ID: 2,29 mm, Wall: 0.86, OD: 4,01 mm**

Art. ID	GE-2.29-PBLK3B-F
Unit	pack of 6 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Tygon R3603/R3607 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Inexpensive all-round tubing for general laboratory applications // Limitations: Potential leaching of plasticizers // Advantages: Transparent, resistant to almost all inorganic chemicals, does not age, low gas permeability, smooth inner surface