

Contour Flared End PVC Tubing - 3 Tag, 95mm, Tag Colour: purple/white, ID: 2,79 mm, Wall: 0.86, OD: 4,51 mm

Art. ID	GE-2.79-PW3B-95-F
Unit	pack of 6 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Tygon R3603/R3607 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Inexpensive all-round tubing for general laboratory applications // Limitations: Potential leaching of plasticizers // Advantages: Transparent, resistant to almost all inorganic chemicals, does not age, low gas permeability, smooth inner surface