

**Common Accessories - Contour Flared End Solva Tubing - 2 Tag, 140mm, Tag Colour: orange/yellow, ID: 0,51 mm, Wall: 0,91, OD: 2,33 mm**

Art. ID	GE-0.51-OY-SF-F
Unit	each
Deliverydetails	No Dangerous Good /not restricted

**Description**

All tubes made from Tygon 4040 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Special tubing for hydrocarbons, petroleum products and distillates // Limitations: Not recommended for strong acids and alkalis, foodstuffs, beverages and medicines. Potential leaching of plasticizers // Advantages: Specially formulated to transport hydrocarbons, petroleum products and distillates, ideal for gasoline, kerosene heating oils, cutting liquids and coolants based on glycols, low gas permeability /// Thermo (VG)|ICP-MS|X Series|Axial