

Common Accessories - Contour Flared End Solva Tubing - 2 Tag, 140mm, Tag Colour: white/white, ID: 1,02 mm, Wall: 0.86, OD: 2,74 mm

Art. ID	GE-1.02-WHT-SF-F
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

Instruments: Thermo (VG) Thermo (TJA) /// All tubes made from Tygon 4040 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Special tubing for hydrocarbons, petroleum products and distillates // Limitations: Not recommended for strong acids and alkalis, foodstuffs, beverages and medicines. Potential leaching of plasticizers // Advantages: Specially formulated to transport hydrocarbons, petroleum products and distillates, ideal for gasoline, kerosene heating oils, cutting liquids and coolants based on glycols, low gas permeability