

**Contour Flared End Solva Tubing - 3 Tag, 95mm, Tag Colour: orange/orange, ID: 0,89 mm, Wall: 0.86, OD: 2,61 mm**

Art. ID	GE-0.89-ORG3B-95-SF-F
Unit	pack of 6 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Tygon 4040 or equivalent // Physical Properties: Thermoplastic, PVC based material with plasticizer // Special Properties: Special tubing for hydrocarbons, petroleum products and distillates // Limitations: Not recommended for strong acids and alkalis, foodstuffs, beverages and medicines. Potential leaching of plasticizers // Advantages: Specially formulated to transport hydrocarbons, petroleum products and distillates, ideal for gasoline, kerosene heating oils, cutting liquids and coolants based on glycols, low gas permeability