

Contour Flared End Tygon MH Tubing - 2 Tag 72mm between tags, Tag Colour: yellow/blue, ID: 1,52 mm, Wall: 0.86, OD: 3,24 mm

Art. ID	GE-1.52-YB-MH3-F
Unit	pack of 6 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Tygon MHSL 2001, Suitable for MIBK // Physical Properties: Special thermoplastic of high purity (without additives, without plasticizer, environmentally-friendly disposal) // Special Properties: Environmentally-friendly tubing especially designed for solvents // Limitations: Limited service-life, cannot be repeatedly sterilized // Advantages: No additives, no plasticizers, suitable for methyl ethyl ketone (MEK), acetone and other corrosive solvents, extremely smooth inner surface, easily disposable