

Contour Flared End Viton Tubing - 2 Tag, 127mm, Tag Colour: orange/orange, ID: 0.89 mm, Wall: 0.85, OD: 2.59 mm

Art. ID	GE-VS0.89-ORG-F
Unit	pack of 6 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Fluran HC F-5500A or equivalent // Physical Properties: fluorpolymerelastomer // Special Properties: Special tubing for concentrated acids and corrosive solvents // Limitations: Limited service-life // Advantages: High chemical resistance, low gas permeability, wide temperature range