

Contour Flared End Viton Tubing - 2 Tag, 127mm, Tag Colour: yellow/yellow, ID: 1.42 mm, Wall: 0.85, OD: 3.12 mm

Art. ID	GE-VS1.42-YLW-F
Unit	pack of 6 tubes
Deliverydetails	No Dangerous Good /not restricted

Description

All tubes made from Fluran HC F-5500A or equivalent // Physical Properties: fluorpolymerelastomer // Special Properties: Special tubing for concentrated acids and corrosive solvents // Limitations: Limited service-life // Advantages: High chemical resistance, low gas permeability, wide temperature range