

**Silicon (Si), 1000 µg/mL in 0,05% HF**

Art. ID AP-CZ90531F5  
Unit 500 mL/HDPE  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	1000	µg/mL		(NH <sub>4</sub> ) <sub>2</sub> SiF <sub>6</sub> : 99,99 %