

Tin (Sn), 1000 µg/mL in 1% HF + 5% HNO₃

Art. ID AP-CZ90551FN5
Unit 500 mL/HDPE
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|-----------------|
| | Tin (Sn) | [7440-31-5] | 1000 | µg/mL | | Sn: 99,999 % |