

3-Amino-1-propanol-d6

Art. ID TRC-A628133-2.5MG
Unit 2.5 mg
Deliverydetails to be checked

Description

Application Note: 3-Amino-1-propanol-d6 is a deuterated version of 3-Amino-1-propanol (A628130) which is used in the synthesis of tetrathiafulvalenes displaying optical and electrochemical properties. Also used in the synthesis of lipid-like materials for the delivery and encapsulation of RNAi therapeutics.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-Amino-1-propanol-d6	[59720-07-9]				