

Unalloyed Copper XI (chips)

Art. ID NIST-454
Unit 35 g
Deliverydetails No Dangerous Good /not restricted

Description

Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	24	µg/g		
	Arsenic (As)	[7440-38-2]	46	µg/g		
	Bismuth (Bi)	[7440-69-9]	19	µg/g		
	Gold (Au)	[7440-57-5]	7,5	µg/g		
	Lead (Pb)	[7439-92-1]	66	µg/g		
	Selenium (Se)	[7782-49-2]	479	µg/g		
	Silver (Ag)	[7440-22-4]	286	µg/g		
	Tellurium (Te)	[13494-80-9]	27	µg/g		
	Tin (Sn)	[7440-31-5]	2,2	µg/g		
	Zinc (Zn)	[7440-66-6]	7	µg/g		
	Copper (Cu)	[7440-50-8]	99,84	%		