

## **Beryllium-Copper (17200) (chips)**

Art. ID	NIST-459
Unit	50 g (chips)
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM®) is in the form of chips sized between 0.50 and 1.18 mm sieve openings (35 and 16 mesh). Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	1,82	%		
	Cobalt (Co)	[7440-48-4]	0,221	%		
	Nickel (Ni)	[7440-02-0]	0,039	%		
	Iron (Fe)	[7439-89-6]	0,079	%		
	Silicon (Si)	[7440-21-3]	0,077	%		
	Aluminium (Al)	[7429-90-5]	0,044	%		
	Tin (Sn)	[7440-31-5]	0,005	%		
	Zinc (Zn)	[7440-66-6]	0,002	%		
	Chromium (Cr)	[7440-47-3]	0,005	%		
	Lead (Pb)	[7439-92-1]	0,001	%		
	Magnesium (Mg)	[7439-95-4]	0,007	%		