

## **Tungsten-Molybdenum EPMA**

Art. ID	NIST-480
Unit	disc
Deliverydetails	No Dangerous Good /not restricted

### Description

This standard is a composite consisting of a core of W-20 wt%Mo alloy wire (1 mm diameter), embedded in pure Mo rod (5 mm diameter) onto which has been electroplated a layer of pure W (1 mm thick). Each standard is 1 mm thick.  
/// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tungsten (W)	[7440-33-7]	78,5	%		
	Molybdenum (Mo)	[7439-98-7]	21,5	%		