

2,4,4'-Trichlorobiphenyl (PCB 28) 50 µg/mL in 2,2,4-Trimethylpentane (Isooctane)

Art. ID AP-CE100-1.5I
Unit 1,5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4,4'-Trichlorobiphenyl I (PCB 28)	[7012-37-5]	50	µg/mL		