

**2,2',5,5'-Tetrachlorobiphenyl (PCB 52) 50 µg/mL in 2,2,4-Trimethylpentane (Isooctane)**

Art. ID AP-CE101-1.5l  
Unit 1,5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',5,5'-Tetrachlorobiphenyl (PCB 52)	[35693-99-3]	50	µg/mL		