

**2,3,3',4,4'-Pentachlorobiphenyl (PCB 105) 50 µg/mL in 2,2,4-Trimethylpentane (Isooctane)**

Art. ID AP-CE120-1.5I  
Unit 1,5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,3',4,4'-Pentachloro biphenyl (PCB 105)	[32598-14-4]	50	µg/mL		