

2,3,3',4',6-Pentachlorobiphenyl (PCB 110) 50 µg/mL in 2,2,4-Trimethylpentane (Isooctane)

Art. ID AP-CE108-1.5I
Unit 1,5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|--------------|---------------------|-------|--------|--------|
| | 2,3,3',4',6-Pentachloro biphenyl (PCB 110) | [38380-03-9] | 50 | µg/mL | | |