

2,3,4,4',5-Pentachlorobiphenyl (PCB 114) 50 µg/mL in 2,2,4-Trimethylpentane (Isooctane)

Art. ID AP-CE121-1.5l
Unit 1,5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,4,4',5-Pentachlorobiphenyl (PCB 114)	[74472-37-0]	50	µg/mL		