

**2,3',4,4',5-Pentachlorobiphenyl (PCB 118) 50 µg/mL in 2,2,4-Trimethylpentane (Isooctane)**

Art. ID AP-CE109-1.5I  
Unit 1,5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3',4,4',5-Pentachloro biphenyl (PCB 118)	[31508-00-6]	50	µg/mL		