

2',3,4,4',5-Pentachlorobiphenyl (PCB 123) 50 µg/mL in 2,2,4-Trimethylpentane (Isooctane)

Art. ID AP-CE122-1.5I
Unit 1,5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2',3,4,4',5-Pentachloro biphenyl (PCB 123)	[65510-44-3]	50	µg/mL		