

2,2',3,3',4,4',5-Heptachlorobiphenyl (PCB 170) 50 µg/mL in 2,2,4-Trimethylpentane (Isooctane)

Art. ID AP-CE112-1.5I
Unit 1,5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,3',4,4',5-Heptachlorobiphenyl (PCB 170)	[35065-30-6]	50	µg/mL		